



Excelarate™ SPD Endo Wipe™

Keep your endoscopes and accessories safe and reprocessing-compliant with our sterile, non-linting wipes.



Excelarate SPD Endo Wipe is AORN®, SGNA®, ASGE®, and AAMI® compliant.

Excelarate Advantages

COMPLIANCE

Safety is a priority in medical equipment reprocessing. Major manufacturing brands such as OLYMPUS®, FUJIFILM®, and PENTAX® all recommend using fresh, non-linting, sterile wipes to dry endoscopes and accessories after disinfection and liquid chemical sterilization.

PATIENT SAFETY

Adherence to Original Equipment Manufacturer (OEM) reprocessing instructions is vital to Patient Safety and the maintenance of your endoscopes. The following are OEM reprocessing instruction excerpts:

- “Sterile materials should be used for all reprocessing steps to avoid recontamination, which may result in infection.” - FUJIFILM
- “After rinsing, thoroughly dry the electrical contacts of the endoscope connector by wiping them with sterile lint-free cloths. Otherwise, hard water residue may be deposited on the electrical contacts, which may cause an abnormal endoscopic image when the endoscope is used.”
- OLYMPUS



DRYING

- Excelarate wipes exhibit outstanding absorption capability. One wipe is enough to dry one endoscope.
- Excelarate wipes are sturdy, making them easier to grip.

OPTIMIZED PACKAGING

- Each easy-open pack contains two folded wipes for compact storage and reduced footprint.
- The sealed packaging protects against dust, contamination, and moisture that may cause bacterial growth.

Ethylene Oxide (EO) Sterilization

- Excelarate Wipes are Ethylene Oxide (EO) sterilized, and each package has an EO indicator for visual confirmation.
- All packets indicate the date of manufacture and sterilization expiration date.

Value

- Excelarate wipes cost approximately 60% less than other sterile non-linting wipes.



OEM USER GUIDELINES

OEM	Point of Use	Manual Cleaning	Drying pre-HLD/LCS	Drying post-HLD/LCS	Accessories
OLYMPUS®	Cloth or Sponges	Cloth,Brushes or Sponges	Cloth	Sterile Cloth	Sterile Cloth
PENTAX®	Gauze Pads	Gauze Pad or Cloth	Gauze Pad or Cloth	Sterile Gauze Pad or Cloth	Sterile Gauze
FUJIFILM®	Gauze	Cloth or Sponge	Gauze	Sterile Gauze	Sterile Gauze

Item Code	Description
SW-001	Excelarate™ SPD Endo Wipe™: Sterilized, Non-Linting Poly Cellulose Wipes, 9" x 9" (23 cm x 23 cm), 400 hundred wipes per case (two wipes per individually sealed pack).

Have a question? Let's Connect.

Customer Service: **877.572.9985**

Facsimile: **877.572.9986**


Email: **orders@clinicalchoice.com**

or visit us on the web at

www.clinicalchoice.com

Clinical Choice LLC

5574 Garden Village Way, Ste. D400
Greensboro, NC 27410



Scan this QR code to
purchase/request a quote for
Excelarate Wipes